### NASA TECH BRIEF



NASA Tech Briefs are issued to summarize specific innovations derived from the U.S. space program, to encourage their commercial application. Copies are available to the public at 15 cents each from the Clearinghouse for Federal Scientific and Technical Information, Springfield, Virginia 22151.

# Spectrophotometric Technique Quantitatively Determines NaMBT Inhibitor in Ethylene Glycol-Water Solutions

#### The problem:

In studies of heat transfer media, problems of corrosion and formation of insoluble residues made it necessary to find an improved method of quantitatively determining the sodium mercaptobenzothiazole (NaMBT) inhibitor in ethylene glycol-water solutions. The improved method must be faster and more accurate than wet chemical techniques. It must also reduce analysis time, require smaller samples, and should be able to detect extremely small concentrations of mercaptobenzothiazole.

#### The solution:

A spectrophotometric method, using a ratio-recording ultraviolet-absorption spectrophotometer, permits analysis of NaMBT in ethylene glycol-water solutions with high accuracy.

#### How it's done:

Analyses are made using a spectrophotometer with a hydrogen lamp energy source. Sample preparation depends upon the condition of the sample and the expected NaMBT concentration; samples which contain insoluble residues or which appear cloudy must be filtered prior to analysis. The samples also must be quantitatively diluted using distilled water (1:100) because of the extreme sensitivity of the technique.

The spectrophotometer records a continuous curve representing the concentration of the NaMBT in the ethylene glycol-water solutions. The samples are scanned from 360 to 200 millimicrons, using a one centimeter sample cell (with distilled water in the reference cell). The maximum absorption at 311 millimicrons is compared to a previously prepared calibration curve to determine the amount of NaMBT present.

This technique has a limit of detection of 250 ppm and is accurate to 0.005 percent.

#### Notes:

- 1. A calibration curve was prepared for use with the spectrophotometer in determining the amount of NaMBT in heat transfer solutions. Several solutions of known NaMBT concentration (in the range of 0.060 to 0.220 percent by weight of 50 percent NaMBT for adequate coverage) were prepared and analyzed. The concentrations were plotted against the absorptance at 311 millimicrons for each sample and the curve was prepared, based on Beer's Law of Absorption.
- Care must be taken to insure that no interferences are present in the ethylene glycol-water solution that will give an absorptance reading at 311 millimicrons.
- Inquiries concerning this innovation may be directed to:

Technology Utilization Officer Manned Spacecraft Center Houston, Texas 77058 Reference: B67-10573

#### Patent status:

No patent ation is contemplated by NASA.

Source: G. G. Garrard of North American Aviation, Inc. under contract to Manned Spacecraft Center (MSC-11496)

Category 03



## MASA TECH BEIER

# Talker to the state of the stat

man and the second of the seco

and the second s

de la digne en recent de la companie de la companie

the land of the contract of the land of th

States August

ASA A winter of the second of the Common of the Common of the control of the cont

Milwood C

Table 1 - The content of a content of the content o

kiton, of all 20 griftspecifies sell cabes to the specifies the cabes being the cabes of February and the specifies the cabes of the cabe of the cabes of the cabes of the cabes of the cabe of